## The solution<br/>from Pender de Balanciafrom Pender de Balanciafrom Pender de Balanciabalancia

## Does your Mammography Information System offer the following solution?

PenAutoReader utilizes adaptive learning and pattern recognition technology to eliminate manual data input for tracking and/or using structured report systems.

With the PenAutoReader system the radiologist continues to transcribe narrative reports while the system automatically extracts exam information from the narrative reports. Many facilities have attached AutoReader to their RIS to save time and labor collecting data.

Pertinent tracking information such as exam type, letter type, recall interval, and BI-RADS are automatically extracted and stored for MQSA audit, tracking, statistics and unresolved case report generation.

The PenAutoReader adapts to each radiologist by each exam type to enhance data extraction. PenAutoReader frees staff and radiologists from entering information for tracking and statistical analysis.

Interested in learning more about PenRad Solutions, call 763.475.3388 for more info or to schedule your personal webinar or visit our website at www.penrad.com

